Caframo

High Torque Stirrer User Manual

Language Index

English.....	1
French.....	17
Spanish.....	33
German	47

Contents

Language Index	1
Online Warranty Registration.....	1
Packing List	2
General Instructions	2
Introduction	3
Safety Instructions	3
Approvals.....	5
Technical Data.....	6
Torque Capabilities.....	7
Dimensions.....	8
Installation	8
Environmental Operating Conditions	9
General Operating Instructions	9
Zero Torque Feature	11
Auto Start	12
Cleaning and Service	12
Warranty and Product Liability	13

Online Warranty Registration

To register the warranty for your Caframo stirrer online, please complete the registration found at

CaframoLabSolutions.com/product-registration/

Packing List


Compare the contents of the shipment to the items on the following list to ensure all parts are received with the stirrer. Do not discard the packing materials until all parts are accounted for. It is recommended that the original packing materials be kept for transportation and storage.

- Overhead Stirrer
- Chuck with Key
- Chuck Guard
- Calibration Certificate
- Instruction manual
- Cord

Cord Sets included with models below:

BDCXXXX

BDCXXXX220U1 (EU,UK)

BDCXXXX220AU

General Instructions

Read Safety Instructions: Unpack the stirrer carefully. Ensure the stirrer is not damaged and verify the contents against the packing list. If the stirrer is damaged or the contents do not correspond with the packing list, notify your supplier immediately.

Read the manual carefully! Make sure that all users read this manual thoroughly before operating the stirrer. For further questions or inquiries, please contact:

Caframo Limited

501273 Grey Road 1
Georgian Bluffs ON N0H 2T0, Canada

Tel: 519-534-1080

Toll Free: 800-567-3556 (within USA and Canada)

Website: CaframoLabSolutions.com

E-mail: labsales@caframo.com

Introduction

All Caframo stirrers are specifically designed for professional use in the laboratory and for use in industry to stir numerous liquids with various properties and applications.

BDC1850 - Rotates at 12-1800 rpm with a maximum torque of 583 Ncm

BDC3030 - Rotates at 20-3000 rpm with a maximum torque of 350 Ncm

BDC6015 - Rotates at 40-6000 rpm with a maximum torque of 175 Ncm

The stirrer uses a high efficiency brushless DC motor. The speed is kept constant with continuous feedback to the motor. All stirrers have integrated overload protection which ensures no damage can occur to the electronics, motor or transmission when an overload situation occurs, such as too great a torque demand or stirrer overheating.

The stirrers incorporate smart logic for automatic transmission control. This logic determines the appropriate speed range at which the stirrer should operate. See Torque Capabilities section. The automatic selection can be overridden by pressing the speed button.

The stirrer is designed for reliability and durability. There are no wear parts that require user service during the normal life of the stirrer.

Safety Instructions

1. Before use, all operators must be familiar with the stirrer and should read this entire manual as labeled on the back of the stirrer. If the stirrer is used with accessories not provided or recommended by the manufacturer, or used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.

-  2. The stirrer must be securely fixed to a stable support. Mount to a Caframo Stand (Part Number A110) with a Caframo clamp (Part Number A120). If other stands or clamps are used, the stirrer must be secure to prevent fall if tipped up to ten degrees from the vertical.

-  3. The stirrer must be supplied with the voltage specified on the serial label.

-  4. **CAUTION:** This is not an explosion proof stirrer. Do not use with highly flammable or explosive materials. Do not operate the stirrer in explosive atmospheres, with hazardous substances or under water.

-  5. Spinning impellers can cause severe personal injuries. Operators must use extreme care and good judgement when mixing at any speed. The stirrers have higher torque capabilities than conventional stirrers.

- ⚠ 6. All mixing paddles and impellers must be in good condition with straight shafts. If the stirrer vibrates at high speeds check the paddle shaft for damage and repair or replace it.
- ⚠ 7. Extreme care must be taken when mixing chemicals so that no chemicals are splashed outside the mixing vessel. Care must be taken when changing to faster mixing speeds. On power up, the stirrer will display its set speed and will ramp to that speed when the start/pause button is pushed. Always start at lowest speed if unsure of maximum safe speed.
- ⚠ 8. Ensure that the mixing impeller does not contact the vessel.
- ⚠ 9. Do not operate while standing in water. Keep the stirrer dry and do not immerse any part, except the mixing paddle, into any liquids.
- ⚠ 10. Ensure that no loose clothing, jewelry, or hair can become entangled in any rotating parts. A fast spinning chuck can cause injury to an operator. Use chuck guard provided when stirrer is in use.
- ⚠ 11. Ensure the power button and mains power cord of the stirrer can be accessed immediately, directly and without risk at any time, and above all, in case of fault.
- ⚠ 12. Shaft rotation can be stopped by pressing the start/pause button or power button. This does not disconnect power to the stirrer's internal drive circuit. In an emergency, power can be interrupted by disconnecting the mains cord.
- ⚠ 13. Wear appropriate personal protective equipment such as eye protection and any clothing suitable for the application.
- ⚠ 14. Repairs must be carried out only by technicians authorized by Caframo.
- ⚠ 15. Stirrer must be used with supplied power cord(s). The cord(s) and grounded plug types are important safety elements. Always connect the stirrer to a grounded outlet. Do not alter the plug.

Approvals

Standards

CAN/CSA C22.2 No. 61010-1-12

CAN/CSA-C22.2 No. 61010-2-051-15

UL 61010-1 (3rd Edition)

UL 61010-2-051

FCC Part 15, Class A

ICES003

IEC 61010-1, IEC 61010-2-051

EMC/EMI standard IEC 61326-1

RoHS	Restriction of Hazardous Substances (EU)
WEEE	Waste Electrical and Electronic Equipment (EU)
CSA	Canadian Standards Association
CE	Communauté Européenne (European Community)
UL	Underwriters Laboratories
FCC	Federal Communications Commission
IEC	The International Electrotechnical Commission
ICES	Interference Causing Equipment Standards
EMC	Electromagnetic Compatibility
EMI	Electromagnetic Interference

1.  Caution/Warning
2.  Recyclable Electric Equipment

Technical Data

Speed Range	BDC1850: 12-1800 rpm BDC3030: 20-3000 rpm BDC6015: 40-6000 rpm
Speed Readout Accuracy	+/- 1 rpm or +/- 1% of reading
Maximum Torque	BDC1850: 580 Ncm (821 in-oz) BDC3030: 350 Ncm (496 in-oz) BDC6015: 175 Ncm (247 in-oz)
Torque Readout Accuracy	+/- 12 Ncm (1 in-lb) or +/- 5% of reading
Voltage	110-240 VAC, 50/60 Hz +/- 10% of the nominal voltage
Installation Category (in accordance with IEC664)	II
Maximum Current Draw	2.6 Amps
Maximum Output Power	1/5 hp, 150 Watts (BLDC Motor)
Weight	5 kg (11 lbs)
Maximum Volume	BDC1850: 80 L (21 Gal) BDC3030: 60L (15.5 Gal) BDC6015: 25 L (6.5 Gal)
Maximum Viscosity	BDC1850: 90,000 cP BDC3030: 50,000 cP BDC6015: 20,000 cP
Housing	Cast aluminum, epoxy coated TENV (totally enclosed non-ventillated)
IP Rating	IP 4X
Keypad	Sealed, chemical resistant
Power Cord Length	150 cm (5 ft) from stirrer to locking connector 50 cm (1.5 ft) from locking connector to plug
Chuck	Fits 13 mm (33/64 in) dia output shaft of stirrer Opens to 10 mm (3/8 in)
Chuck Guard	Length 67 mm, (2.64 in) Diameter 54 mm, (2.13 in) Molded clear chemical resistant polyurethane
Support Rod	12.7 mm (0.5 in) dia 140 mm (5.25 in) long
Drive	2 speed transmission
Through Shaft	Accepts mixing shafts up to 10 mm dia (3/8 in)


Torque Capabilities

Stirrer	Low Speed				High Speed			
	Min Speed rpm	Max Speed rpm	Peak Torque Ncm	Cont. Torque* Ncm	Min Speed rpm	Max Speed rpm	Peak Torque Ncm	Cont. Torque* Ncm
1850	12	360	583	320	117	1800	120	75
3030	20	600	350	190	195	3000	72	43
6015	40	1200	175	85	389	6000	36	21


*Continuous torque at 50% of rpm range and 25 degrees Celcius.

Torque Conversion Table		
In-lb to In-oz	Multiply by 16	50 in-lb = 800 in-oz = 568 Ncm
In-lb to Ncm	Multiply by 11.36	30 in-lb = 400 in-oz = 341 Ncm
Ncm to In-oz	Multiply by 1.416	15 in-lb = 240 in-oz = 170 Ncm


Low Speed

High Speed

Dimensions

Installation

The stirrer must be securely fixed to a stable support - e.g. clamped to an original Caframo stand that is on a stable, level surface. See Safety Instructions.

Chuck removal and attachment

Align chuck set screw with pilot hole on the motor shaft. Slide chuck onto motor shaft. Tighten the set screw using the hex end of the provided chuck key.

Impeller attachment

Open the chuck fully using the provided chuck key. Insert the impeller shaft into the open chuck and adjust the impeller to the desired height. The hollow shaft allows the impeller shaft to come through the top of the stirrer if necessary. Tighten the chuck's outer ring by hand and then secure with the chuck key using all three insertion points.

Environmental Operating Conditions

The stirrer must operate in the following conditions:

Indoors

1. Altitudes up to 2000 m (6500 ft).
2. Temperatures from 5°C to 40°C (41°F to 104°F).
3. Maximum relative humidity 80% for temperatures up to 31°C (88°F) decreasing linearly to 50% relative humidity at 40°C (104°F).
4. Pollution degree 2 in accordance with IEC 664.


General Operating Instructions

-  Make sure that the stirrer is securely mounted as specified in the installation instructions. Also make sure that the impeller is securely attached in the desired position.
-  Connect the power cord to a grounded supply.
-  The stirrer is now ready for use to mix liquids and liquid/solid solutions. The operator, ensuring that all safety instructions are followed must determine appropriate impeller and container specifications.
-  The stirrer can operate continuously if the torque demand falls below the continuous torque operation line (see Torque Capabilities). When using the torque curve charts, should the ambient temperature be above 25°C, then the 40°C continuous operation line must be used.
-  If the stirrer is used for an extended period of time above the appropriate continuous temperature operation line (see Torque Capabilities), then it may shut down. The stirrer may be started after it has cooled down.

Button Operating Functions

-  Power - Turns the stirrer on and off. When power button is pressed for the first time stirrer is activated in “pause” mode. The set speed is displayed. When the Power button is pressed a second time, the stirrer is turned off. The last set speed is remembered.
-  Start/Pause - stops the stirrer when operating. Start/Pause returns the stirrer to the previous condition and previous display or starts the stirrer if stopped.
-  Display - When running, toggles display between speed (rpm) and torque (Ncm).
-  Speed - determines speed range (high/low speed) and overrides automatic choice when allowable. The LEDs indicate which speed range is operating.
-  Increase/Decrease rpm Buttons - Press buttons once to increase / decrease set speed by one rpm. Press and hold to increase / decrease speed at a greater rate. As the stirrer is auto ramping to a faster speed the user can stop the increase in speed by pressing the decrease button. This changes the set speed to the present actual value.

Keypad Membrane

Indicator Lights

- High Speed / Low Speed LEDs - These indicate the stirrer's set range.
- Overload LED - This indicator will flash when an overload condition is approached. This indicator will remain on when an overload condition is reached and the stirrer is shut down.

- Power LED is on when stirrer is on.
- Pause LED is on when stirrer is paused and set speed is displayed.
- Display Button - The LEDs indicate which parameter is displayed. When paused, the display will always show set speed. The rpm LED will be flashing indicating this is a set speed.

Overload conditions will occur when:

- The torque demand is outside the stirrer's capabilities (solution is too viscous). Overload LED will illuminate. Reduce torque requirements of the application to continue mixing.
- The stirrer overheats due to extended operation at very high torque and/or speed demands. The stirrer completely shuts down and must cool before operation can resume.

Zero Torque Feature

Theory

In mixing applications the magnitude of torque depends on three parameters: First, the speed in revolutions per minute; second, the diameter of the propeller; and third, the viscosity of the material being mixed. Additionally, a small amount of torque is required to move the mixing impeller. In small bench top applications this additional or overhead torque is negligible and is generally independent of the force required to mix material.

Zero Torque Feature

The zero torque feature allows the operator to null out the additional overhead torque of the stirring apparatus similar to a tare feature on a scale. The feature can also be used during the mixing process and allows the operator to see change in torque caused by the change in material viscosity.

Run the stirrer at low speed for a half hour before zeroing torque. As the motor warms up the torque can change slightly. Setting the torque when the stirrer is warm will increase the accuracy.

A negative number indicates a reduction in viscosity and a positive number incates an increase in viscosity.

Other button functionality

- If the stirrer is powered off when in zero torque mode (Ncm LED flashing), the stirrer will be in this mode when powered back on. If the stirrer is unplugged or power is interrupted, the stirrer is reset and the feature will not be active.
- The speed can be changed any time - while the stirrer is running or in pause. When the speed is changed, the speed in rpm is displayed while the speed is adjusted. If torque was displayed prior to the change in speed, the display will

change back to displaying torque once the new speed is set.

- The zero torque feature can only be activated when the stirrer is running.
- When the zero torque button is pressed while the stirrer is running and rpm is displayed, actual torque will be displayed. Press the zero torque button again, the Ncm LED flashes indicating the stirrer is in zero torque mode.

Auto Start

This stirrer includes an optional Auto Start feature. When set, the stirrer will resume mixing after a power failure, when power has been removed by control equipment (i.e., PLC) or when the AC power has been cut (i.e., stirrer unplugged). When power is restored, the stirrer will automatically ramp up to the last set speed and resume mixing. The stirrer ships from the factory with Auto Start disabled.

To set the Auto Start feature first ensure the stirrer is powered off. Press and hold the Increase Button then press the Power Button. When active, Auto Start is indicated by a solid LED dot in the bottom right corner of the display.

To disable Auto Start first ensure the stirrer is powered off. Press and hold the Decrease Button then press the Power Button.

Cleaning and Service

Cleaning

Use a mild soap and water solution and wipe the housing and control panel of the stirrer with a damp cloth. Do not use chlorine bleach, chlorine-based cleanser, abrasives, ammonia, steel wool, or scouring pads with metal content or similar harsh solvents or abrasives. These may damage the surface of the stirrer.

Service

Do not attempt to service or repair stirrer. If the stirrer housing is opened, the warranty becomes void. Contact Supplier for return authorization and return instructions. Clean the stirrer so that it is free from any materials which may be harmful. Provide an SDS or other safety data sheets for all chemicals used in conjunction with this stirrer where appropriate. Place the stirrer and its parts into the original packing to protect and prevent damage during transport. Follow shipping directions provided by Caframo when return is authorized and RMA issued.

CAUTION:

 Failure to clean, service, and handle the stirrer as outlined above can cause it to fail which could impair the safety of the user.

Warranty and Product Liability

Warranty

Caframo warrants the herein described product of its manufacture for three years from receipt of the equipment by the purchaser, against defects in material and workmanship. Warranty registration is available on the Service and Support page of our website CaframoLabSolutions.com.

This limited warranty covers parts and labour. In the event of a warranty claim, contact Caframo. If the cause is determined to be a manufacturing fault, Caframo will repair or replace all defective parts to restore the stirrer to operation.

This warranty may only be altered by a specifically published amendment. No individual has authorization to alter the provisions of this warranty policy or its amendments. Expendable items are not covered by this warranty. Damages due to corrosion, accidental breakage or submersion are not covered.

Product Liability

Under no circumstances shall Caframo be liable for indirect, consequential or special damages of any kind. Caframo reserves the right to make technical changes without prior notice.

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IEC 61010-1, IEC 61010-2-051

EMC/EMI standard IEC 61326-1

RoHS	Restriction of Hazardous Substances (EU)
WEEE	Waste Electrical and Electronic Equipment (EU)
CSA	Canadian Standards Association
CE	Communité European (European Community)
UL	Underwriters Laboratories
FCC	Federal Communications Commission
IEC	The International Electrotechnical Commission
ICES	Interference Causing Equipment Standards
EMC	Electromagnetic Compatibility
EMI	Electromagnetic Interference

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