

Torque You Trust

Robust 568 Ncm torque
for precision in intense
mixing tasks

Precision

Toggle display to view
RPM or Ncm to track
changes in viscosity

Easy Adjustments

Through Shaft for quick
impeller and batch
changes

No Downtime

Continuous viscous
mixing without fail
3 year warranty

1800 RPM

Dispersion

Pigment Grind
High Shear Mixing

Formulation

Let-Down
Low Shear
Blending Prop

Adhesive Systems

1K or 2K, Easy-Clean
PTFE Anchor

**12 RPM
568 Ncm**

Canadian
Made



The Kit:

	Ultra Torque Overhead Stirrer	Safety Stand	Fastening Clamp	Impeller Shaft	Dispersion Blade	Pitched Blade Propeller	PTFE Anchor
Caframo Part	BDC1850	A110	A120	A722	A164	A531	A183